



USB DIGITAL MICROSCOPE CAMERA

- This excellent microscope consists of a colour zoom CMOS video camera mounted in a microscope styled housing.
- Connects to a PC via USB connection.
- Up to 640 x 480 image resolution.
- 5x optical zoom, with maximum magnification x 130
- 30fps maximum video frame rate.

CODE	DESCRIPTION
QC-3244	USB Digital Microscope Camera System



STEREO/DISSECTING MICROSCOPE

- High quality unit with all the practical features.
- Dual magnification (2X & 4X) with 10X wide field adjustable eye pieces.
- Built in illumination for both transmitted & incident light.
- 240 Volt operation.
- Extension for viewing large items.
- Plastic cover.
- Clear, frosted & black stage options.
- Clear image viewing.

A MICROSCOPE THAT ALL FISH FARMERS COULD USE

CODE	DESCRIPTION
ST-3C	Stereo/Dissecting Microscope



MICROSCOPE INCLINED

- Built in halogen lamp 6v-20W and solid state transformer.
- Superior double condenser.
- Mechanical stage.
- Objectives: Achromatic 10x, 40/45x (sl), 100x (sl) oil immersion.
- Eye piece: WF 10x and H 5x or 10x or 15x
- Comes with polishing cloth, dust cover and working instructions in thermacole box.

CODE	DESCRIPTION
XSP-102A	Microscope inclined



MULTI-FUNCTION MICROSCOPE

- Multi-function biological microscope.
- Binocular professional microscope
- Objectives 4X, 10X, 40X & 100X (oil) objectives plus 10X & 16X eye pieces.
- Built in adjustable light intensity.
- Double layer adjustable mechanical stage with measurements.
- Sliding binocular head & focus for easy viewing.
- Spare halogen globe, dust cover & immersion oil.

IDEAL FOR THE SERIOUS AQUACULTURIST, SCHOOL, UNIVERSITY OR COLLEGES, AND RESEARCH INSTITUTES

CODE	DESCRIPTION
BF202	Multi-Function Bio-Microscope

MICROSCOPE COVER SLIPS

- 22mm x 22mm plain square glass slips.
- Suitable for use with the MS50 slides.



CODE	DESCRIPTION
MCS100	Microscope Cover Slips pkt of 100