

## MICROSCOPE SLIDES

- Clear glass slides with ground edges.
- Slides are 25mm x 75mm x 1.2mm thick.
- Pack of 50.



CODE	DESCRIPTION
MS50	Microscope Slides pkt of 50

## MICROSCOPE CASE

- Easy storage of slides
- Ideal protection for your microscope slides.
- Slides and cover slips not included.

CODE	DESCRIPTION
MS90	Microscope Case



## DISSECTING KIT

- Two dissecting needles.
- Probe with eye.
- Scissors sharp/sharp.
- Scissors sharp/blunt.
- Forceps blunt/sharp straight.
- Forceps sharp - curved.
- Scalpel plus handle with interchangeable blade.
- Hook with chain.

*"This will last you for years"*

CODE	DESCRIPTION
DK	Dissecting kit in Zipped Case



## DISSECTING TRAY

- Inexpensive, polyethylene dissecting trays.
- Includes a long-lasting vinyl pad.
- Great for use in labs and with juvenile fish.

CODE	DESCRIPTION
DT9003	Dissecting Tray



## MEASURING CYLINDERS

- Moulded polypropylene transparent flasks.
- Autoclavable and chemically resistant.
- Permanently moulded graduations.

CODE	DESCRIPTION
MC50	Measuring Cylinder 50mL
MC100	Measuring Cylinder 100mL
MC250	Measuring Cylinder 250mL



## ERLENMEYER FLASKS

- Moulded polypropylene transparent flasks.
- Autoclavable and chemically resistant.

CODE	DESCRIPTION
EF250	Erlenmeyer Flask 250mL
EF500	Erlenmeyer Flask 500mL

