

SYRINGES

CODE	DESCRIPTION
SY1	1ml Syringe - 0.1m Graduations
SY25	25ml Syringe - 1ml Graduations

TRANSFER PIPETTES

- Made from polypropylene
- 0.5ml moulded graduations for easy use.

CODE	DESCRIPTION
TP20	3ml Graduated Transfer Pipettes pkt of 20
TP50	3ml Graduated Transfer Pipettes pkt of 50

PETRI DISHES

These disposable Petri dishes are 90mm x 15mm in size and made from polystyrene.

CODE	DESCRIPTION
PD20	Petri Dishes pkt of 20

LATEX GLOVES

- Strong latex gloves in large size.

CODE	DESCRIPTION
LATEX	Latex Disposable Gloves (100/pkt, Large)

WEIGHING PANS

CODE	DESCRIPTION
WP2	Medium Weighing Pan (10 pack)
WP3	Large Weighing Pan (10 pack)

EZY POURER 90ML DISPENSER

Easy to use, screw the Ezy Pourer onto your chemical bottle and squeeze until the Ezy Pourer fills to the selected level and pour. When liquid enters the Ezy Pourer and hits the adjustable level indicator filling ceases so that the exact amount of liquid is left in the chamber for you to pour out, its that simple.

CODE	DESCRIPTION
EP90	Ezy Pourer 90ml Dispensing Bottle Top

BACTERIA FILTER

- Reusable or throw away, filters air for use with single species algae culture, artemia and rotifers, bio-assay work, experiments, etc.
- Filtering down to 0.3 Microns, 6-8mm Hose barb connection

Code	Description
FL2	0.3micron Bacteria Air Filter

