

PLANKTON NETS, MINIATURE

- Plankton nets have a 12.7cm diameter net mouth and measure 38cm long.
- Net mouth has stainless steel bridle ring with 3 point tow line attachments.
- Nets include a 125mL sampling bottle.
- Great for shallow water sampling.



LM007

CODE	DIAMETER (cm)	LENGTH (cm)	MESH-SIZE (microns)
LM007	12.7	38	153
LM008	12.7	38	80

PLANKTON NETS

- Designed for biological plankton & oceanographic studies.
- Used for a variety of aquaculture applications.
- Includes a stainless steel mouth ring/towing bridle, cod end ring.
- Rubber coated hose clamp.
- 1 Litre capacity cod end jar.
- Effective conical design with a large mouth allows for easy collecting of large samples.



PKN1

CODE	DIAMETER (cm)	LENGTH (cm)	MESH SIZE (microns)
PKN1	30	90	35
PKN2	30	90	63
PKN3	30	90	80
PKN4	30	90	100
PKN5	30	90	150
PKN6	50	150	80
PKN7	50	150	150
PKN8	50	150	200
PKN9	50	150	335
PKN10	50	150	500
PKN11	75	225	250

